

LCAM startup manual, Andor microscope (A2.42)

update: Mar 2023, RB



Contact person: Ronald Breedijk, C2.267 (+ 7860), R.M.P.Breedijk@uva.nl


2nd contact person: Mark Hink, C2.264 (+ 6211), M.A.Hink@uva.nl

Information microscope: <http://www.lcam-fnwi.nl/facilities/andor>



*Before using the Andor microscope, the user (and co-workers) should have had the official intake discussion with the LCAM-staff, succeeded the LCAM-confocal training course & exam and had an individual training at the Andor microscope. **Bookings** can be requested via cam.microscopy@gmail.com*

Start-up procedure Andor confocal

1. Switch on the powersocket at the optical table (between lasers and camera) using the red switch.
2. Switch on the computer: user: **imcb** password: **Microscopy@23**
3. If needed switch on the Oko Touch temperature controller and CO₂ regulator.



CO₂ switch (on the left side of remote control) and temperature switch (on the backside)





4. If needed switch on the Nikon Intensilight fluorescence lamp.

5. Start software, NIS Elements (from desktop) and Select OK for ""Andor""




6. Start Laser control software (from desktop)

7. Switch on all the laser-lines you will need in the laser program from 7a->7d



7c. Select laser line(s)
7d. Select Run button




7a. Select Card: 3

7b. Press shutter button

8. Data storage

Never store data on the local harddisk, all data should be stored at the dataserer. Data present on the local PC will be deleted without further notice. Be aware that the storage of data on the sever will be your own responsibility as well. Although there is a regular backup of the server we will not take any responsibility for lost or damaged data, so make backups yourself. Contact Mark Hink in order to get access to the data-server from your office computer.


Login using the **green Connect dataserver** icon at the desktop: Type your userID (UvAnetID), group and password. After succesful login a network drive U:\ will be visible where you should store your data.



Basic handling of the Andor microscope

Put a droplet of the correct immersion liquid on top of the objective (air, water or oil). A small bottle of **Nikon** oil can be found nearby. In case the oil bottle is empty: **NEVER use the immersion oil from other brands (Leica/Zeiss/Olympus)** but contact Ronald or Mark for a refill.

1. Open the **NIS Elements** window.
2. For the control of the microscope the **Ti Pad** window can be used. If it is not open already you can open it by a right mouse click having the mouse cursor positioned at the grey background of the NIS Elements main screen. Then select **Acquisition controls** and then **Ti Pad**.
3. Select the optical configuration you would like to use.
4. Switch from port 1 to port 0 and back to port 1 (software bug).
5. Optimize Laserpower, exposure time for optimal image quality.



Ti Pad

Nosepiece
100x 100x 60x 20x 10x 60x

Escape
Escape Z

Light Path
E100 On Memory Recall
L100 R100 L80 Focus Offset Dichroic Mirror In

Z Drive
Move by step[μm]: 0.1 1 10 20.0
Z[μm]: 3742.2
Accuracy[μm]: 0.000

Lamps
DIA 3.0 12.0 3.0 [V]
TIRF 0.0 7000.0 0.0
Coarse Fine Extra Fine

Filters
Turret2 Turret1
Excitation Wheel

Condenser
Zoom 1.00x

Configure...

AOTF Pad

405 nm 440 nm 488 nm 514 nm
561 nm 593 nm 640 nm

1: 405 nm 0 [%]
2: 440 nm 0 [%]
3: 488 nm 20 [%]
4: 514 nm 0 [%]
5: 561 nm 0 [%]
6: 593 nm 0 [%]
7: 640 nm 0 [%]

Stimulation
405 nm Power 100 [%]
Manual Test Shot
Pulse Test Shot Pulse time: 100 [msec]

Shutters
AOTF Configure...

Filters, Shutters and Switchers
CSU AOTF Aux1
Turret2 Turret1
CSU DM
Excitation Wheel
WF

DU-897 1 Settings

Format For Live No Binning
Format For Capture No Binning
Auto Exposure 100 ms
Readout Mode EM Gain 10 MHz at 14-bit
EM Gain Multiplier 100
Conversion Gain 5.1x
Temperature -70.8 °C *Desired temp. differs!*

Commands

DU-897 2 Settings

Format For Live No Binning
Format For Capture No Binning
Auto Exposure 200 ms
Readout Mode EM Gain 10 MHz at 14-bit
EM Gain Multiplier 50
Conversion Gain 5.1x
Temperature -69.5 °C

Commands

CSU Pad

CSU
1500 5000 5000 [rpm] Sync
CSU DM

ND Acquisition

Experiment: ND Acquisition DU-897 2

Save to File Record Data...
Custom Metadata

Order of Experiment Timing...
Time schedule

Phase	Interval	Duration	Loops
#1	10 sec	5 min	31

Close Active Shutter when Idle Perform Time Measurement (0 ROIs)
Switch Transmitted Illuminator off when Idle (1.00 s)

Autofocus None At the Beginning Define...
Execute before each Time
Execute after each Time
Advanced for: Time Phase 1 Apply to All
Execute before Time Phase
Execute after Time Phase

XYZ Navigation

XY [μm] Z1 [μm]
Coarse Fine
XY step 0.500 FOV 100%
Z2 [μm]
0.1 1 10 5.000



XY - Ti XYDrive
X [μm]: 1531.100 Range: <-55.000, 55.000> mm Move
Y [μm]: 1032.100 Range: <-37.500, 37.500> mm
Z1 - Ti ZDrive
Z [μm]: 3742.23 Range: <0.0, 10000.0> μm Move
Escape Refocus

Z2 - NIDAQ Piezo Z (name: Piezo Z)
Z [μm]: 50.000 Range: <0.000, 100.000> μm Move
Piezo

XY=[1531.100, 1032.100]μm,
Z=3692.23μm, Z1=3742.23μm, Z2=50.00μm

LUTs
315 G: 1.00 1843

6. Setting up a Z-stack in the **Capture Z-series** window (see right figure) .



7. Setting up a timelapse sequence in the **Capture Timelapse** menu (See left figure)

Components of Andor microscope

- Objectives:**
 - 20x air ELWD NA 0.45
 - 40x oil NA 1.3
 - 60x TIRF oil NA 1.49
 - 100x TIRF oil NA 1.49
- Filters for visual inspection of fluorescence -> see optical schematic
- Lasers:** 405 nm, 445 nm, 488 nm, 514 nm, 561 nm, 593 nm and 640 nm

Extra features Andor microscope

1. The Andor allows to use a **spinning disk**. This can be used for confocal images. The advantage over standard confocal microscopy is the imaging speed can be much faster.
2. The Andor can also be used for TIRF imaging such that you only “illuminate” the border part of a cell, till ca 100 nm above the glass.



Switching off procedure of the Andor microscope

1. Shut down NIS software
2. Logoff using the **red Disconnect dataserver** icon at the desktop
3. Switch off powersocket
4. Switch off Hg lamp
5. Switch off Okolab environmental control (if used)